

MEDUTEK



KATASPIR 20

Respiratory tract suction



KATASPIR 20

Respiratory tract suction
efficient & reliable

KATASPIR 20 is suitable for nasal, oral and endotracheal aspiration of bodily fluids from adults and children. Indicated for tracheotomized patients, minor surgical applications and post-operative therapy. KATASPIR 20 finds application in outpatient and inpatient care, elderly care and in private care.

Efficient and **oilless and maintenance-free** piston pump. Also intended for non-stop operation.

It can be used with FLOVAC disposable liners or reusable containers.

SCOPE OF DELIVERY

- ✓ Small and handy unit
- ✓ Proven technology
- ✓ AquaSafe filter
- ✓ MEDUTEK service-package
- ✓ Silicone tubing set
- ✓ Reusable jar 1000 ml with cover

MAIN APPLICATIONS

ENT, Home Care, Room/Ward, Tracheostomy

Art.-no.: 483-KAT20NM





TECHNICAL SPECIFICATIONS



- ✓ Max. free air flow rate: 16 l/min.
- ✓ Max. vakuüm (*adjustable*): - 0,75 bar / - 75 kPa
- ✓ Low noise level: 60,5 dBA



- ✓ 2 years of warranty
- ✓ 5 years MEDUTEK service-package



SCOPE OF DELIVERY:

- ✓ KATASPIR 20
- ✓ Silicon tubing set with filter
- ✓ 1000 ml jar with cover
- ✓ Instruction manual



- ✓ Duty cycle: Non-stop operation
- ✓ Motor: Oilless and maintenance-free piston pump
- ✓ Power feeding: 230V-50Hz



- ✓ Size: 350 x 210 x 180 mm
- ✓ Weight: 2,5 KG



Art.-no.: 483-KAT20NM

 +49-(0)421 427 480

 INFO@MEDUTEK.COM



SUCTION UNITS



SUCTION UNITS

 +49-(0)421 427 480

 INFO@MEDUTEK.COM