



KATASPIR 20 eco

The compact suction unit



KATASPIR 20 eco

The compact suction unit
small and handy

KATASPIR 20 eco is suitable for nasal, oral and endotracheal aspiration of bodily fluids from adults and children. Indicated for tracheotomized patients, minor surgical applications and post-operative therapy. KATASPIR 20 eco finds application in outpatient and inpatient care, elderly care and in private care.

Efficient and **oilless and maintenance-free** piston pump. Also intended for non-stop operation.

It can be used with FLOVAC disposable liners or reusable containers.

SCOPE OF DELIVERY

- ✓ Small and handy unit
- ✓ Proven technology
- ✓ AquaSafe filter
- ✓ MEDUTEK service-package
- ✓ Silicone tubing set
- ✓ Reusable jar 1000 ml with cover

MAIN APPLICATIONS

ENT, Home Care, Room/Ward, Tracheostomy





TECHNICAL SPECIFICATIONS



- ✓ Max. free air flow rate: 15 l/min.
- ✓ Max. vakuüm (*adjustable*): - 0,75 bar / - 75 kPa
- ✓ Low noise level: 59,6 dBA



- ✓ 2 years of warranty
- ✓ 5 years MEDUTEK service-package



SCOPE OF DELIVERY:

- ✓ KATASPIR 20 eco
- ✓ Silicon tubing set with filter
- ✓ 1000 ml jar with cover
- ✓ Instruction manual



- ✓ Duty cycle: Non-stop operation
- ✓ Motor: Oilless and maintenance-free piston pump
- ✓ Power feeding: 230V-50Hz



- ✓ Size: 250 x 190 x 160 mm
- ✓ Weight: 2,2 KG



Art.-no.: 483-KAT20ECO

+49-(0)421 427 480


INFO@MEDUTEK.COM



SUCTION UNITS



SUCTION UNITS

 +49-(0)421 427 480

 INFO@MEDUTEK.COM