

MEDUTEK



KATASPIR 30

The all-round suction pump



KATASPIR 30

The all-round suction pump
power & high-performance

KATASPIR 30 is suitable for nasal, oral and endotracheal aspiration of bodily fluids (mucus or catarrh) from adults and children. Indicated for tracheotomized patients, surgical applications and post-operative therapy, KATASPIR 30 finds application in out-patient and in-patient care, elderly care, in private care and at professional level due to its powerful suction capacity of 40 l/min.

Efficient and **oilless and maintenance-free** piston pump. Also intended for non-stop operation.

It can be used with FLOVAC disposable liners or reusable containers.

SCOPE OF DELIVERY

- ✓ Small and handy unit
- ✓ Proven technology
- ✓ AquaSafe filter
- ✓ MEDUTEK service-package
- ✓ Silicone tubing set
- ✓ Reusable jar 1000 ml with cover

MAIN APPLICATIONS

Endoscopy, Gastroenterology, General Surgery, Home Care, Obstetrics, Room/Ward, Tracheostomy, Veterinary





TECHNICAL SPECIFICATIONS



- ✓ Max. free air flow rate: 40 l/min.
- ✓ Max. vakuüm (*adjustable*): - 0,80 bar / - 80 kPa
- ✓ Low noise level: 60,5 dBA



- ✓ 2 years of warranty



SCOPE OF DELIVERY:

- ✓ KATASPIR 30
- ✓ Silicon tubing set with filter
- ✓ 1000 ml jar with cover
- ✓ Instruction manual



- ✓ Duty cycle: Non-stop operation
- ✓ Motor: Oilless and maintenance-free piston pump
- ✓ Power feeding: 230V-50Hz



- ✓ Size: 350 x 210 x 180 mm
- ✓ Weight: 3,6 KG



Art.-no.: 483-KAT30

+49-(0)421 427 480

INFO@MEDUTEK.COM



SUCTION UNITS



SUCTION UNITS



+49-(0)421 427 480



INFO@MEDUTEK.COM