

MEDUTEK



KATASPIR 36BR (Battery)

Mobility with high-power



KATASPIR 36BR

Mobility with high-power

KATASPIR 36BR is a portable suction unit for the aspiration of body liquids, oral, nasal and tracheal aspiration in adults or children. For professional use and with a powerful aspiration of max 36 l/min, KATASPIR 36BR is lightweight for easy transport. Three different options for operation: AC/DC adapter, internal rechargeable battery and 12V make it a versatile suction unit approved for NON-STOP operation without overheating. Provided with visual and acoustic alarm indicating low battery level.

Efficient and **oilless and maintenance-free** piston pump. Also intended for non-stop operation.

It can be used with FLOVAC disposable liners or reusable containers.

SCOPE OF DELIVERY

- ✓ Small and handy unit
- ✓ Proven technology
- ✓ AquaSafe filter
- ✓ MEDUTEK service-package
- ✓ Silicone tubing set
- ✓ Reusable jar 1000 ml with cover

MAIN APPLICATIONS

Ambulance Car, Emergency, Home Care, Patient Transfer, Room/Ward, Tracheostomy, Veterinary

Art.-no.: 483-KAT36BR

 +49-(0)421 427 480

 INFO@MEDUTEK.COM



SUCTION UNITS



TECHNICAL SPECIFICATIONS



- ✓ Max. free air flow rate: 36 l/min.
- ✓ Max. vakuüm (*adjustable*): - 0,80 bar / - 80 kPa
- ✓ Low noise level: 65,5 dBA



- ✓ 2 years of warranty



SCOPE OF DELIVERY:

- ✓ KATASPIR 36BR
- ✓ Silicon tubing set with filter
- ✓ 1000 ml jar with cover
- ✓ Instruction manual
- ✓ Power Supply Unit (AC/DC adapter, 12V adapter)



- ✓ Duty cycle: Non-stop operation (by mains only)
- ✓ Motor: Oilless and maintenance-free piston pump
- ✓ Power feeding: 14V 4A by AC/DC Universal Adapter, Internal Rechargeable Pb Battery, - Autonomy 60min - Recharging Time 240 min - Car Adapter



- ✓ Size: 350 x 210 x 180 mm
- ✓ Weight: 4,39 KG

Art.-no.: 483-KAT36BR

+49-(0)421 427 480

INFO@MEDUTEK.COM



SUCTION UNITS



SUCTION UNITS

 +49-(0)421 427 480

 INFO@MEDUTEK.COM