

*A company from Bremen with more than 40 years of experience
in healthcare!*

CATALOGUE BREATHING TUBES



BREATHING

CPAP-Breathing Tubes

Professional Equipment for CPAP, Sleep Apnea, Respiratory Therapy and Homecare

Who we are

The Bremen based company MEDUTEK Handelsgesellschaft für Medizintechnik mbH & Co. KG is regarded as a reliable partner by our customers and suppliers of medical technology for over 40 years.

MEDUTEK has many years of experience in the field of medical devices and their accessories and is committed to the individual interests and wants for its worldwide partners. The company maintains a certified quality management system according to EN ISO 13485:2016.

Breathing tubes in different variety

MEDUTEK offers a wide range of high-quality and certified breathing tubes:

- Our tubes come in variable diameters, lengths, colors and tube material
- Our breathing tubes are cost-effective and economically advantageous
- Immediate availability and prompt delivery of all of our tube styles
- Customized tubes can be ordered
- All our tubes are CE certified and meet the ISO standards

Our hoses are also used during surgery for the effective extraction of surgical smoke.

Order process:

On the following pages we will give you an overview of our breathing tubes.

In addition, we can also supply customised breathing tubes on request. Please do not hesitate to contact us.

Our item numbers are created as follows:

Our product code-style number-inner diameter-length-colour

For example: Item no. 181-603191800GR

181: Our Product code

603: Style 603 Tube

19mm: Diameter of the tube

1800 mm: Length of the tube

GR: Grey color

Grey: GR; Transparent: CL; Black: SCH; Beige;








Style 303
Autoclavable



Style 305
Autoclavable

MEDUTEK Style 303/305 tubing is lightweight, flexible, autoclavable, and a very good alternative to silicone tubing. The Style 305 tube offer integrated cuffs with the same inner diameter as the tube.

Category	Breathing Tube Style 303	Breathing Tube Style 305
Applications	Breathing Circuit, Respiratory Care, CPAP, Home Care, Ventilator	Breathing Circuit, Respiratory Care, CPAP, Home Care, Ventilator, Fume extraction during surgery
Other Features	<ul style="list-style-type: none"> - sterilizable - reusable - smooth interior - high stability and flexibility - custom length 	<ul style="list-style-type: none"> - sterilizable - reusable - smooth interior - high stability and flexibility - custom length
Inner Diameter	19 mm	10 mm bis 32 mm
Length Range	150 – 3000 mm	150 – 3000 mm
Tube Material	TPC-ET, Hytrel	TPC-ET, Hytrel
Cuffs	integrated	integrated
Available Color		 
Max. Temperature	-20°C to 140 °C	-20°C to 140 °C
Compatibility	Cuffs with ISO-22 mm connector	Cuffs with the same inner diameter
Regulatory	CE, ISO 10993, ISO 5367, REACH-RoHS compliant	CE, ISO 10993, ISO 5367, REACH-RoHS compliant
Item-Nr.	181-303XXXXXXXX	181-305XXXXXXXX

CONTACT
US FOR
QUANTITY
DISCOUNT!





Style 403
Autoclavable



Style 405
Autoclavable

MEDUTEK Style 403/405 tubing is intended for use with ventilators and respiratory therapy devices for adults. Style 403/405 breathing tubes are supplied with cuffs designed for 22 mm ISO conical connectors. Style 403 and 405 differ in their inner diameter.

Category	Breathing Tube Style 403	Breathing Tube Style 405
Applications	Breathing Circuit, Reusable, Respiratory Care, CPAP, Home Care, Ventilator	Breathing Circuit, Reusable, Respiratory Care, CPAP, Home Care, Ventilator
Other Features	<ul style="list-style-type: none"> - autoclavable - reusable - easy cleaning - smooth interior - high stability and flexibility - custom length 	<ul style="list-style-type: none"> - autoclavable - reusable - easy cleaning - smooth interior - high stability and flexibility - custom length
Inner Diameter	19 mm	10 mm, 15 mm
Length Range	150 – 3000 mm	150 – 3000 mm
Tubing Material	TPC-ET, Hytrel	TPC-ET, Hytrel
Cuffs	Silicone	Silicone
Available Color		
Max. Temperature	-20°C to 140 °C	-20°C to 140 °C
Compatibility	Cuffs with ISO-22 mm connector	Cuffs to Fit
Regulatory	CE, ISO 10993, ISO 5367, REACH-RoHS compliant	CE, ISO 10993, ISO 5367, REACH-RoHS compliant
Item-Nr.	181-403XXXXXXXX	181-405XXXXXXXX



**WHOLESALE
PRICE
ONLY!**



Style 503

Standard



Style 505

Standard

Medutek Style 503/505 Breathing tubes are made of flexible, ribbed, biocompatible, medical-grade plastic material and are available in variable lengths from 150- 3000 mm. The tubes are supplied with cuffs designed for 22 mm ISO conical connectors. Style 503 is available only in 19 mm inner diameter and Style 505 with variable inner diameter.

Category	Breathing Tube Style 503	Breathing Tube Style 505
Applications	Breathing Circuit, Reusable, Respiratory Care, CPAP, Home Care, Ventilator	Breathing Circuit, Reusable, Respiratory Care, CPAP, Home Care, Ventilator
Other Features	<ul style="list-style-type: none"> - easy cleaning - smooth interior - high stability and flexibility - custom length 	<ul style="list-style-type: none"> - easy cleaning - smooth interior - high stability and flexibility - custom length
Inner Diameter	19 mm	10 mm bis 51 mm
Length Range	150 – 3000 mm	150 – 3000 mm
Tubing Material	Polyolefin plastomer and low-density polyethylene	Polyolefin plastomer and low-density polyethylene
Cuffs	integrated	integrated
Available Colors		
Max. Temperatur	-20°C to 70 °C	-20°C to 70 °C
Compatibility	Cuffs with ISO-22 mm connector	Custom
Regulatory	CE, ISO 10993, ISO 5367, REACH-RoHS compliant	CE, ISO 10993, ISO 5367, REACH-RoHS compliant
Item-Nr.	181-503XXXXXXXX	181-505XXXXXXXX

**MEDUSoft
OUR
BEST
SELLING
TUBES!**



Style 603
Standard



Style 605 /
MEDUSoft

MEDUTEK Style 603/605 breathing tubes are supplied with cuffs designed for 22 mm ISO conical connectors. Style 603 and 605 differ in their inner diameter. Style 605 is very well-known in the German market as MEDUSoft.

Category	Breathing Tube Style 603	Breathing Tube Style MEDUSoft
Applications	Breathing Circuit, Reusable, Respiratory Care, CPAP, Home Care, Ventilator	Breathing Circuit, Reusable, Respiratory Care, CPAP, Home Care, Ventilator
Other Features	<ul style="list-style-type: none"> - easy cleaning - smooth interior - high stability and flexibility - custom length 	<ul style="list-style-type: none"> - easy cleaning - smooth interior - high stability and flexibility - custom length
Inner Diameter	19 mm	15 mm
Length Range	150 – 3000 mm	150 – 3000 mm
Tubing Material	Polyolefin plastomer and low-density polyethylene	Polyolefin plastomer and low-density polyethylene
Cuffs	TPE	TPE
Available Color	○ ●	●
Max. Temperatur	-20°C bis 70 °C	-20°C bis 70 °C
Compatibility	Cuffs with ISO-22 mm connector	Cuffs with ISO-22 mm connector
Regulatory	CE, ISO 10993, ISO 5367, REACH-RoHS compliant	CE, ISO 10993, ISO 5367, REACH-RoHS compliant
Item-Nr.	181-603XXXXXXXX	181-605XXXXXXXX



Style 623
Standard



Style 625
Standard


The MEDUTEK non-heated CPAP tube Style 623/625 delivers a continuous positive airway pressure to ensure that the patient is continuously receiving constant flow of pressurised air during sleep. The pressure remains at a constant level during inhalation and exhalation and can be adjusted according to device type and patient's needs. Style 623 and 625 differ with their inner diameter and length.

Category
Applications
Other Features
Inner Diameter
Length Range
Tubing Material
Cuffs
Available Colour
Max. Temperature
Compatibility
Regulatory
Item-Nr.

Breathing Tube Style 623
Respiratory Care, CPAP, Home Care
<ul style="list-style-type: none"> - low flow resistance - high pressure resistance - easy cleaning - smooth inner surfaces - high stability and flexibility
19 mm
1800 mm
PE/PP
PE/PP
○ ●
≤ 60 °C
conical cuffs with ISO 5356-1 connector
CE, ISO 5356
061-62318CL 061-62318GR

Breathing Tube Style 625
Respiratory Care, CPAP, Home Care
<ul style="list-style-type: none"> - low flow resistance - high pressure resistance - easy cleaning - smooth inner surfaces - high stability and flexibility
15 mm
1800 mm
PE/PP
PE/PP
●
≤ 60 °C
conical cuffs with ISO 5356-1 connector
CE, ISO 5356
061-625151800GR

FOR
MORE
INFORMATION
CALL US
NOW!




Hybernite

The Hybernite Rainout Control System is used to prevent condensation of humidified air delivered to patients requiring positive airway pressure or mask ventilation. Hybernite can be used with almost all CPAP devices.

Customized Tubing

Customized tubes from Medutek are the solution to your specific needs. According to your requirement, we can deliver tubings and cuffs made of TPC-ET, LDPE Helix, EVA, Silicone etc. Also the size, length and diameter and shape of the tubes and cuffs can be adjusted.

Category	
Applications	Breathing Circuit, Reusable, Respiratory Care, CPAP, Home Care, Ventilator, Low Pressure Applications
Other Features	<ul style="list-style-type: none"> - heatable tube - easy cleaning - cuffs with humidifier connection - smooth interior - high stability and flexibility - ergonomic hand grips
Inner Diameter	19 mm
Length Range	1830 mm
Tubing Material	PP/PE
Cuffs	PP/PE
Available Colours	
Max. Air Temperatur	5°C to 30 °C
Power Supply	Frequency: 47–63 Hz; 100 V–240 V
Compatibility	Cuffs with ISO-22 mm connector
Regulatory	CE, ISO 8185:2007, ISO 5367:2014
Item-Nr.	061-183019

Accessories

- Tubing connectors
- Bacterial filter
- Fullface mask with or without silicone cushion
- Nasal mask with or without silicone cushion
- Hybernite Power cord
- Hybernite Power supply adapter

Please feel free to contact us if you need any further information.

www.medutek.com



INFO@MEDUTEK.COM



+49-(0)421 427 480



For orders or product information, please contact our Sales Support directly:

Medutek Handelsgesellschaft für Medizintechnik mbH & Co. KG
Stresemannstraße 54
D-28207 Bremen

✉ sales@medutek.com
☎ +49 (0)421 - 427 48-0

MEDUTEK Product Category



Scan the QR Code to reach out for us:



You can find more accessories and alternative devices here:
Scan code with a QR code scanner. Open link.

Other Services that we offer for our partners:

- Product training
- Personal consultation
- Videos in Youtube
- Information material like product data sheet, brochures etc.



MEDUTEK

Handelsgesellschaft für Medizintechnik mbH & Co. KG
Stresemannstraße 54
D-28207 Bremen

Tel.: +49 (0)421- 427 48-0

Email: info@medutek.com

www.medutek.com

ND-MEDUTEK©2022

All rights reserved, errors and omissions excepted.

All content and images in association with this catalogue are owned by MEDUTEK. The information in this catalogue has been provided by MEDUTEK for general information purposes only, we make no representations or warranties of any kind about the accuracy, reliability and suitability of the information.



Follow us on:

